

SCHWARZ BIOSCIENCES	FORMULATION & TECHNOLOGY REPORT	Date
		Project No. Incontinence
TITLE TDS for the treatment of incontinence, part III: Delivery of SPM8224, the free base of Fesoterodine CONFIDENTIAL		Page 1 of 9
		Report No. 30

1. SUMMARY

The report describes in vitro skin permeation characteristics of transdermal delivery systems (TDS) containing SPM8224, the free base of Fesoterodine. Test samples were prepared by either lab-scale solvent coating or hot-melt processing. Patches were tested by means of flux rates across hairless mouse skin, selected samples were also investigated in the LACDR human skin model.

High flux rates of SPM8224 across mouse skin were obtained for all batches. The highest initial flux was achieved with a silicone based hot melt formulation.

Permeation across human skin demonstrated the promising potential of SPM8224, the free base of Fesoterodine, for the treatment of overactive bladder. Based on these in vitro data patches with sizes in the range of 15 to 30 cm² could theoretically delivery 4 to 8 mg/24 h which is the current range of the oral Fesoterodine formulation. The data have to be confirmed in vivo.

Besides this, new acrylic based hot melt adhesive from National Starch & Chemical were evaluated. One very promising formulation could be identified. Provided that other success criteria, such as physico-chemical compatibility with the hot melt process and with different drugs, are met these formulations might fill the gap in hot melt pressure sensitive adhesives for transdermal systems.

Distribution: Original PH DOK F&T, PHA, TS, TT, PH REG, IPM (AS) Summary only/PCD, PH TOX, BA, MOBI, SIL, ILF			
Key words: Fesoterodine, SPM8224, skin permeation in vitro, mouse skin, human skin, hot melt acrylics			
Author	Name	Signature	Date
Head of TS	Dr. A. Breitenbach	<i>A. Breitenbach</i>	
Reviewed by			
Head of TT	Dr. H.-M. Wolff	<i>H.-M. Wolff</i>	
Approved by			
Head of F&T	M.C.F. Hannay	<i>M.C.F. Hannay</i>	

SCHWARZ//BIOSCIENCES	FORMULATION & TECHNOLOGY REPORT	Date
		Project No. Incontinence
TITLE TDS for the treatment of incontinence, part III: Delivery of SPM8224, the free base of Fesoterodine CONFIDENTIAL	Page 2	of 9
	Report No.	30

CONTENTS

	PAGE
1. SUMMARY	1
2. INTRODUCTION AND OBJECTIVES	3
3. MATERIALS AND METHODS	4
4. RESULTS AND DISCUSSION	5
5. CONCLUSION	9


APPENDIX A (Certificates of analysis)

APPENDIX B ()

SCHWARZ BIOSCIENCES	FORMULATION & TECHNOLOGY REPORT	Date
		Project No. Incontinence
TITLE TDS for the treatment of incontinence, part III: Delivery of SPM8224, the free base of Fesoterodine CONFIDENTIAL		Page 3 of 9
		Report No. 30

2. INTRODUCTION AND OBJECTIVES

The objective of the study was to investigate the feasibility of transdermal delivery of SPM8224, the free base of Fesoterodine (scheme 1).



Scheme 1: Chemical structure of the free base of Fesoterodine

Therefore, several lab scale patch batches containing SPM8224 were prepared and investigated by means of in vitro drug permeation across hairless mouse skin. Subsequently, selected samples were investigated in a human skin model, as well.

SCHWARZ BIOSCIENCES	FORMULATION & TECHNOLOGY REPORT	Date
		Project No. Incontinence
TITLE TDS for the treatment of incontinence, part III: Delivery of SPM8224, the free base of Fesoterodine CONFIDENTIAL		Page 4 of 9
		Report No. 30

3. MATERIALS AND METHODS

For a detailed description of the experiments refer to the batch documentation.

Hot-melt patches (exemplary): 8 g of a preformed silicone adhesive were weighed into a beaker and tempered at 160°C for ca. 20 min to achieve a homogenous melt. 0.5 g of inner phase polymer (e.g. poly(ethylene oxide) and 1.5 g of drug were added. After tempering at 160°C for additional 5 min the mixture was homogenized manually and further processed on the pre-tempered Chill-Roll (120°C, 250 µm) for lamination.

5 cm² patches were isolated by manual punching followed by determination of the average patch weight (n=10). Finally, patches were sealed individually in pouches.

Mouse Skin Model (PHA): according to OBU0469.ABV100, rev. 00 (1998) with an active diffusion area of 2.55 cm², a phosphate buffer acceptor phase at pH 6.2 and a temperature of 32°C, n=3

Human Skin Model (LACDR):

according to H. Tanojo et al. , J. Control Rel. 45 (1997) 41-47.

skin from abdomen with a thickness of approx. 250 µm, flux experiment: acceptor phase: PBS, pH= 6.2, temperature: 32°C, diffusion cells with spiral groove (8 cells), groove area: 0.552 cm², dialysis membrane used as separator between skin and acceptor phase
flux: 5 ml/hour PBS, experiment runs for 72 hours, sampling cycle: 3 hours

Analytical Methods (PHA): refer to certificates of analysis

Data Analysis: sigmoidal Boltzmann and linear fit: Microcal Origin 6.0

SCHWARZ/BIOSCIENCES	FORMULATION & TECHNOLOGY REPORT	Date
		Project No. Incontinence
		Page 5 of 9
TITLE TDS for the treatment of incontinence, part III: Delivery of SPM8224, the free base of Fesoterodine CONFIDENTIAL		Report No. 30

4. RESULTS AND DISCUSSION

The free base of Fesoterodine is an oily substance and was used with a purity of approx. 90%. Several lab scale patch batches with different compositions were prepared by either solvent coating or hot melt processing and the preparation presented no difficulties. The properties of these batches are summarized in table 1.

Table 1: Properties of the patch batches

No	Lot No. (Ch.B.)	PSA	Theo. drug loading [% (w/w)]	Patch weight (n=10) [g/m²]
1	20111080	SC acrylic	15	100
2	20111085	HM EVA	15	84
3	20111086	HM silicone	15	63
4	20111087	HM SxS	15	89
5	20111095	HM acrylic 01	15	73
6	20201027	HM acrylic 02	15	121
7	20201028	HM acrylic 03	15	115

SC Acrylic = solvent coating, acrylic type adhesive, Duro Tek 387-2287, National Starch & Chemical (NSC)
 HM EVA = hot melt, ethylene vinyl acetate co-polymer adhesive, Depolix 213, NSC
 HM Silicone = hot melt silicone based adhesive, BioPSA + 5% (w/w) Ozokerite wax, DowCorning
 HM SxS = hot melt, styrene block co-polymer, in house formulation
 HM Acrylic 0x = hot melt acrylic type adhesives, experimental formulations from NSC, refer to Annex

All batches were tested in the hairless mouse skin model.

SCHWARZ BIOSCIENCES	FORMULATION & TECHNOLOGY REPORT	Date
		Project No. Incontinence
		Page of 6 9
TITLE TDS for the treatment of incontinence, part III: Delivery of SPM8224, the free base of Fesoterodine CONFIDENTIAL		Report No. 30

Fig. 1 outlines the cumulative permeation of SPM8224 across hairless mouse skin (HMS) from patches prepared with the most common adhesive classes, acrylic based, ethylene vinyl acetate based (EVA), silicone based and poly(styrene) based (SxS). The flux rates were expressed as permeation of active metabolite, which is the hydrolysis di-hydroxy product (SPM7605).




Fig. 1: Cumulative drug permeation (calculated as permeation of active metabolite, SPM7605)

In all cases high flux rates were observed following a non-linear release kinetics. The flux rates from EVA and SxS matrices were nearly identical. Higher initial rates were obtained for the silicone and the acrylic based systems. While the silicone type batch showed the highest drug permeation in the first 6 to 8 h, the subsequent higher values for the acrylic based patch were mainly due to a higher patch weight.

SCHWARZ/BIOSCIENCES	FORMULATION & TECHNOLOGY REPORT	Date
		Project No. Incontinence
		Page 7 of 9 Report No. 30
TITLE TDS for the treatment of incontinence, part III: Delivery of SPM8224, the free base of Fesoterodine CONFIDENTIAL		

Fig. 2 illustrates the drug permeation across hairless mouse skin (HMS) from new types of pressure sensitive adhesives (PSA): for the first time it was possible to test the performance of acrylic based patch batches which were prepared by hot melt processing. The three different experimental PSA were exclusively obtained from National Starch & Chemical.




Fig. 2: Cumulative drug permeation (calculated as permeation of active metabolite, SPM7605)

While two of these new PSA (batches 20111095 and 20201028) yielded lower flux rates, it was possible with the third PSA (batch 20201027) to achieve flux rates comparable to the solvent born system (batch 20111080). Although still in an experimental stage this new class of acrylic based PSA seems to be capable of closing the final gap in hot melt processing, since up to now the use of acrylics was limited to solvent coating, only.

SCHWARZ//BIOSCIENCES	FORMULATION & TECHNOLOGY REPORT	Date
		Project No. Incontinence
		Page 8 of 9 Report No. 30
TITLE TDS for the treatment of Incontinence, part III: Delivery of SPM8224, the free base of Fesoterodine CONFIDENTIAL		

Two batches were investigated in the LACDR human skin model. Since the skin supporting silicone membrane in this model was hindering drug flux it was exchanged with a dialysis membrane. Fig. 3 demonstrates the obtained flux rates across the composite of excised human skin supported with a dialysis membrane (EHS/DM).




Fig. 3: Cumulative drug permeation (calculated as permeation of active metabolite, SPM7605)

After short lag-times of 6 to 8 h (which are not predictive for the in vivo situation) both batches showed a high steady state flux for at least 2.5-d. From these results patches with a size of already 20 cm² could theoretically deliver therapeutic doses of the free base of Fesoterodine, compared to the oral formulation.

SCHWARZ BIOSCIENCES	FORMULATION & TECHNOLOGY REPORT	Date
TITLE TDS for the treatment of incontinence, part III: Delivery of SPM8224, the free base of Fesoterodine CONFIDENTIAL		Project No. Incontinence
	Page 9 of 9 Report No. 30	

CONCLUSIONS

Several lab scale patch batches containing the free base of Fesoterodine were prepared and tested by means of in vitro drug permeation across mouse skin and human skin.

Based on the results obtained, the flux rates were found to be sufficient for the treatment of overactive bladder with patch sizes in the range of 15 to 30 cm² (equal to ca. 4 to 8 mg/24 h).

Moreover, in cooperation with National Starch & Chemical a new series of hot melt acrylics could be developed. These new adhesives could close an actual gap in hot melt pressure sensitive adhesive (PSA) formulation for transdermal systems. To further evaluate this class of PSA an investigation on physico-chemical compatibility with the hot melt process and different drug substances has to be performed.

ANNEX 1

Copies of the Certificates of Analysis

(signed originals stored at PH DOK)

Analysenzertifikat in vitro Freisetzung durch Mäusehaut

Präparat : SPM 907 Ch.-B.: 20111085
 Wirkstoff : SPM 8224
 Sollgehalt : 15 % TDS - Fläche: 5 cm²
 ABV vom : analog OB 0469.ABV.100 Analysendatum :

Ausgangsanalyse

Bemerkungen: Puffer pH 6,2 KT Freisetzungsfläche 2,55 cm²
 Mäusehautdicke: 1=174 µm; 2=159µm, 3=190µm, 35,5g
 Alter lebend: 8 Wochen, TK-Schrank 10 Wochen SKH-1 ♂
HPLC-Bedingungen:
 Eluent: 600 VT Wasser : 400 VT ACN : 1 VT Trifluoressigsäure
 Temperatur: 35°
 Wellenlänge: 220 nm
 Fluss: 1 ml
 Säule: Waters Spherisorb Nitrite 5µm

Tabelle der kumulierten Freisetzung in µg SPM 8224/1 cm²

Zeit [h]	1	2	3	MW	SD
2	58,74	72,97	104,26	78,66	23,3
4	121,08	130,31	162,74	138,04	21,9
6	188,37	183,76	239,35	203,82	30,8
8	238,21	223,26	286,98	249,48	33,3
24	410,84	379,08	463,27	417,73	42,5
30	470,84	428,55	518,49	472,63	45,0
48	572,15	507,03	629,89	569,69	61,5

MW = Mittelwert SD = Standardabweichung

Achsenabschnitt (b) = 92,7 µg
 Regressionskoeffizient (m) = 14,19 µg/h
 Korrelationskoeffizient (r) = 0,9641 3-24h

Tabelle der kumulierten Freisetzung in µg SPM 7605/1 cm²

Zeit[h]	1	2	3	MW	SD
2	14,65	18,20	20,52	17,79	3,0
4	32,93	34,65	34,64	34,07	1,0
6	52,29	50,29	53,65	52,07	1,7
8	67,36	52,25	66,03	61,88	8,4
24	163,98	136,16	146,39	148,84	14,1
30	187,17	156,71	164,01	169,30	15,9
48	246,26	207,59	207,87	220,57	22,2

MW = Mittelwert SD = Standardabweichung

Achsenabschnitt (b) = 12,1 µg
 Regressionskoeffizient (m) = 5,77 µg/h
 Korrelationskoeffizient (r) = 0,9961 3-24h

Datum

Fabert (PHA)

Dr. Seiffert (PHA)

Analysenzertifikat in vitro Freisetzung durch Mäusehaut

Präparat :	SPM 907 TDS	Ch.-B.:	20111086
Wirkstoff :	SPM 8224		
Sollgehalt :	15%	TDS - Fläche:	5 cm ²
ABV vom :	analog OB 0469.ABV.100	Analysendatum :	

Ausgangsanalyse

Bemerkungen: Puffer pH 6,2 KT Freisetzungsfäche 2,55 cm²
 Mäusehautdicke: 1=188 µm; 2=172µm, 3=167µm, 32,7g
 Alter lebend: 9 Wochen, TK-Schrank 7 Wochen SKH-1 ♂
HPLC-Bedingungen:
 Eluent: 600 VT Wasser; 400 VT ACN; 1 VT Trifluoressigsäure
 Temperatur: 35°
 Wellenlänge: 220 nm
 Fluss: 1 ml
 Säule: Waters Spherisorb Nitrile 5µm

Tabelle der kumulierten Freisetzung in µg SPM 8224/1 cm²

Zeit [h]	1	2	3	MW	SD
2	151,6	86,7	72,8	103,7	42,0
4	264,1	194,6	165,7	208,1	50,6
6	329,6	274,1	234,0	279,2	48,0
8	370,9	340,3	291,7	334,3	39,9
24	438,3	538,7	479,3	485,5	50,5
30	442,8	556,4	503,5	500,9	56,8
48	444,3	563,8	521,3	509,8	60,6

MW = Mittelwert SD = Standardabweichung

Achsenabschnitt (b) =	150,3	µg
Regressionskoeffizient (m) =	14,98	µg/h
Korrelationskoeffizient (r) =	0,9227	3-24h

Tabelle der kumulierten Freisetzung in µg SPM 7605/1 cm²

Zeit[h]	1	2	3	MW	SD
2	28,0	18,2	18,3	21,5	5,6
4	52,4	37,3	38,3	42,7	8,5
6	69,4	51,7	53,3	58,1	9,8
8	82,4	64,3	66,5	71,1	9,8
24	127,3	122,2	138,2	129,2	8,1
30	131,6	130,6	149,8	137,3	10,8
48	135,5	138,8	162,2	145,5	14,6

MW = Mittelwert SD = Standardabweichung

Achsenabschnitt (b) =	24,8	µg
Regressionskoeffizient (m) =	4,51	µg/h
Korrelationskoeffizient (r) =	0,9763	3-24h

Analysenzertifikat

in vitro Freisetzung durch Mäusehaut

Präparat : SPM 907 Ch.-B.: 20111087
 Wirkstoff : SPM 8224
 Sollgehalt : 15 % TDS - Fläche: 5 cm²
 ABV vom : analog OB 0469.ABV.100 Analysendatum :
Ausgangsanalyse

Bemerkungen: Puffer pH 6,2 KT Freisetzungsfläche 2,55 cm²
 Mäusehautdicke: 1=170 µm; 2=162µm, 3=160µm, 31,3g
 Alter lebend: 8 Wochen, TK-Schrank 10 Wochen SKH-1 ♂
HPLC-Bedingungen:
 Eluent: 600 VT Wasser : 400 VT ACN : 1 VT Trifluoressigsäure
 Temperatur: 35°
 Wellenlänge: 220 nm
 Fluss: 1 ml
 Säule: Waters Spherisorb Nitrile 5µm

Tabelle der kumulierten Freisetzung in µg SPM 8224/1 cm²

Zeit [h]	1	2	3	MW	SD
2	44,72	59,11	64,84	56,22	10,4
4	85,65	140,79	132,93	119,79	29,8
6	121,63	218,35	205,37	181,78	52,5
8	150,10	269,83	254,75	224,89	65,2
24	257,35	417,04	435,07	369,82	97,8
30	317,34	466,51	490,29	424,71	93,7
48	384,12	541,07	603,89	509,69	113,2

MW = Mittelwert SD = Standardabweichung

Achsenabschnitt (b) = 77,2 µg
 Regressionskoeffizient (m) = 12,87 µg/h
 Korrelationskoeffizient (r) = 0,9521 3-24h

Tabelle der kumulierten Freisetzung in µg SPM 7605/1 cm²

Zeit[h]	1	2	3	MW	SD
2	17,57	20,90	16,42	18,30	2,3
4	33,57	45,48	32,25	37,10	7,3
6	49,28	70,25	50,44	56,66	11,8
8	62,80	89,55	64,09	72,15	15,1
24	137,51	190,84	138,64	155,66	30,5
30	167,04	214,82	156,88	176,25	33,4
48	204,06	268,98	203,75	225,60	37,6

MW = Mittelwert SD = Standardabweichung

Achsenabschnitt (b) = 15,3 µg
 Regressionskoeffizient (m) = 5,98 µg/h
 Korrelationskoeffizient (r) = 0,9908 3-24h

Datum

Fabert (PHA)

Dr.Seiffert (PHA)

Analysenzertifikat in vitro Freisetzung durch Mäusehaut

Präparat :	SPM 907 TDS	Ch.-B.:	20111095
Wirkstoff :	SPM 8224		
Sollgehalt :	15%	TDS - Fläche:	5 cm ²
ABV vom :	analog OB 0469.ABV.100	Analysendatum :	

Ausgangsanalyse

Bemerkungen: Puffer pH 6,2 KT Freisetzungsfläche 2,55 cm²
 Mäusehautdicke: 1=182 µm; 2=169 µm, 3=176 µm, 32,9 g
 Alter lebend: 9 Wochen, TK-Schrank 7 Wochen SKH-1 ♂
HPLC-Bedingungen:
 Eluent: 600 VT Wasser : 400 VT ACN : 1 VT Trifluoressigsäure
 Temperatur: 35°
 Wellenlänge: 220 nm
 Fluss: 1 ml
 Säule: Waters Spherisorb Nitrite 5µm

Tabelle der kumulierten Freisetzung in µg SPM 8224/1 cm²

Zeit [h]	1	2	3	MW	SD
2	25,7	18,1	21,6	21,8	3,8
4	59,7	47,3	51,8	52,9	6,3
6	87,4	69,2	77,4	78,0	9,1
8	110,2	94,2	100,3	101,6	8,1
24	199,6	212,8	211,4	207,9	7,3
30	253,2	245,3	245,5	248,0	4,5
48	322,0	313,3	318,8	318,1	4,4

MW = Mittelwert SD = Standardabweichung

Achsenabschnitt (b) =	22,4	µg
Regressionskoeffizient (m) =	7,96	µg/h
Korrelationskoeffizient (r) =	0,9849	3-24h

Tabelle der kumulierten Freisetzung in µg SPM 7605/1 cm²

Zeit [h]	1	2	3	MW	SD
2	11,9	10,9	10,7	11,2	0,6
4	26,9	26,7	24,3	26,0	1,4
6	40,0	38,9	35,6	38,2	2,3
8	51,3	52,1	45,7	49,7	3,5
24	112,7	131,2	121,5	121,8	9,3
30	151,7	155,3	140,9	149,3	7,5
48	198,4	210,1	189,2	199,2	10,5

MW = Mittelwert SD = Standardabweichung

Achsenabschnitt (b) =	6,5	µg
Regressionskoeffizient (m) =	4,87	µg/h
Korrelationskoeffizient (r) =	0,9965	3-24h

Analysezertifikat in vitro Freisetzung durch Mäusehaut

Präparat :	SPM 907 TDS	Ch.-B.:	20201027
Wirkstoff :	SPM 8224		
Sollgehalt :	15%	TDS - Fläche:	5 cm ²
ABV vom :	analog OB 0469.ABV.100	Analysendatum :	

Ausgangsanalyse:

Bemerkungen: Puffer pH 6,2 KT Freisetzungsfäche 2,55 cm²
Mäusehautdicke: 1=150 µm; 2=148 µm, 3=158 µm, 29,4g
Alter lebend: 8 Wochen, TK-Schrank 12 Wochen SKH-1 ♂
HPLC-Bedingungen:
Eluent: 600 VT Wasser : 400 VT ACN : 1 VT Trifluoressigsäure
Temperatur: 35°
Wellenlänge: 220 nm
Fluss: 1 ml
Säule: Waters Sphertsorb Nitrile 5µm

Tabelle der kumulierten Freisetzung in µg SPM 8224/1 cm²

Zeit [h]	1	2	3	MW	SD
2	51,3	109,4	80,4	80,4	29,0
4	129,0	217,5	232,3	192,9	55,9
6	208,8	300,8	308,0	272,5	55,3
8	261,7	364,5	422,5	349,6	81,5
24	549,2	591,4	966,0	702,2	229,4
30	634,2	633,6	1081,0	782,9	258,1
48	801,5	680,0	1236,8	906,1	292,8

MW = Mittelwert SD = Standardabweichung

Achsenabschnitt (b) = 87,0 µg
Regressionskoeffizient (m) = 26,42 µg/h
Korrelationskoeffizient (r) = 0,9836 3-24h

Tabelle der kumulierten Freisetzung in µg SPM 7605/1 cm²

Zeit [h]	1	2	3	MW	SD
2	14,1	20,5	17,3	17,3	3,2
4	32,2	41,9	56,9	43,7	12,5
6	49,3	59,5	73,1	60,6	12,0
8	65,2	77,4	93,1	78,6	14,0
24	180,8	217,9	213,9	204,2	20,3
30	216,1	248,6	256,0	240,2	21,2
48	284,6	290,7	313,5	296,3	15,2

MW = Mittelwert SD = Standardabweichung

Achsenabschnitt (b) = 8,4 µg
Regressionskoeffizient (m) = 8,24 µg/h
Korrelationskoeffizient (r) = 0,9978 3-24h

Analysenzertifikat

in vitro Freisetzung durch Mäusehaut

Präparat :	SPM 907	Ch.-B.:	20201028
Wirkstoff :	SPM 8224		
Sollgehalt :	15%	TDS - Fläche:	5 cm ²
ABV vom :	analog OB 0469.ABV.100	Analysendatum :	

Ausgangsanalyse

Bemerkungen: Puffer pH 6,2 KT Freisetzungsfäche 2,55 cm²
 Mäusehautdicke: 1=158 µm; 2=168µm, 3=171µm, 31,4g
 Alter lebend: 8 Wochen, TK-Schrank 12 Wochen SKH-1 ♂
HPLC-Bedingungen:
 Eluent: 600 VT Wasser ; 400 VT ACN : 1 VT Trifluoressigsäure
 Temperatur: 35°
 Wellenlänge: 220 nm
 Fluss: 1 ml
 Säule: Waters Spherisorb Nitrile 5µm

Tabelle der kumulierten Freisetzung in µg SPM 8224/1 cm²

Zeit [h]	1	2	3	MW	SD
2	67,0	46,4	9,0	40,8	29,4
4	166,9	120,0	23,2	103,4	73,3
6	240,8	177,4	33,6	150,6	106,2
8	294,4	217,3	39,8	183,8	130,6
24	314,8	282,5	41,1	212,8	149,6
30	370,7	318,0	43,0	243,9	176,0
48	445,8	418,3	48,1	304,1	222,1

MW = Mittelwert SD = Standardabweichung

Achsenabschnitt (b) =	84,4	µg
Regressionskoeffizient (m) =	6,13	µg/h
Korrelationskoeffizient (r) =	0,7917	3-24h

Tabelle der kumulierten Freisetzung in µg SPM 7605/1 cm²

Zeit[h]	1	2	3	MW	SD
2	18,4	13,9	6,7	13,0	5,9
4	40,8	33,6	15,1	29,9	13,3
6	60,8	51,1	21,3	44,4	20,8
8	80,3	67,0	26,9	58,1	27,8
24	282,3	256,8	86,4	208,5	106,5
30	347,9	294,7	93,9	245,5	134,0
48	444,3	372,3	117,1	311,2	171,9

MW = Mittelwert SD = Standardabweichung

Achsenabschnitt (b) =	-8,1	µg
Regressionskoeffizient (m) =	8,96	µg/h
Korrelationskoeffizient (r) =	0,9990	3-24h

Diffusion experiment with SPM 907 patches on full human skin

Experiment number:

907DF004

Purpose of the study:

To investigate the permeation rate of SPM 907 through skin from two newly formulated patches.

Patch:

Active ingredient: SPM 907

Batch numbers: 20111085 and 20111087

Patch area: 5 cm²

Active ingredient content: 15% SPM 8224

Skin donor:

Birth date: 1994

Sex: female

Skin from: abdomen (belly)

Thickness of dermatomised skin: approximately 250 µm

Diffusion experiment:

Date:

Used cells:

diffusion cells with spiral groove (n=8); groove area: 0.552 cm²

Separator between acceptor phase and skin/patch:

Diachema dialysis membrane, type 10.14, supplied by Diener, München, Germany.

Manufactured from neutral cellulose, molar weight cut-off: 5000; thickness (dry): 25 µm.

Pretreated according to the manufacturer's recommendations.

Diameter of separator, skin and patch punch-cores: 1.8 cm

Setup diffusion cells:

Cell nr.	Cells with/without skin	Batch
1 & 2	with	20111085
3 & 4	with	20111087
5 & 6	without	20111085
7 & 8	without	20111087

Acceptor phase:

PBS pH=8.2

Measured temperature waterbath:

32.0

°C

Flux of acceptor phase:

5

ml/hour

Total duration of the experiment: 72 hours, samples are collected in 3 hour periods.

Observations during dermatomisation, cell assembly, disassembly, etc.

Some of the glue from the patch disks destined for cells 7 and 8 remained on the protective foil.

Mass and volume data on the collected fractions

measured density of the used acceptor phase:

1,007	g/ml
-------	------

1,812

Faktor zur Umrechnung auf cm³:

Flux time (hours)	cell nr.	mass (tubes (g)) empty	full	volume fractions (ml)	8224 g/ml	Fraction 1,812	Milligramm	DIOH g/ml	Fraction 1,812	g/ml	Fraction 1,812	Milligramm
3	1	16,948	32,389	16,329	0,00	0,000	0,000	0,02	0,327	0,592		
	2	17,037	32,958	15,810	0,00	0,000	0,000	0,01	0,159	0,286		0,439
	3	17,103	32,928	15,715	0,00	0,000	0,000	0,02	0,314	0,570		
	4	17,148	32,959	15,703	0,00	0,000	0,000	0,02	0,314	0,569		0,569
	5	17,010	33,141	15,703	0,00	0,000	0,000	0,02	0,314	0,569		0,569
	6	17,147	32,885	15,419	11,90	160,626	345,412	0,84	13,456	24,362		27,688
	7	17,058	32,362	15,410	14,83	200,073	416,893	1,11	17,105	30,995		27,688
	8	17,132	32,816	15,217	19,96	303,589	550,103	1,51	22,978	41,637		45,638
6	1	16,632	33,032	16,268	0,01	0,163	0,235	0,10	1,627	2,947		
	2	16,782	32,711	16,518	0,00	0,000	0,000	0,147	0,01	0,158	0,287	1,617
	3	16,689	32,654	15,755	0,01	0,158	0,235	0,09	1,418	2,569		
	4	17,020	32,867	15,737	0,01	0,157	0,235	0,19	2,560	5,416		3,994
	5	17,187	33,320	16,021	8,16	130,730	238,853	0,60	8,613	17,418		18,077
	6	17,126	32,667	15,433	5,93	81,518	185,830	0,87	10,340	18,796		
	7	16,897	32,273	15,269	6,88	105,062	180,353	0,54	6,245	14,941		
	8	17,104	32,806	15,563	8,18	127,550	231,120	0,83	9,823	17,800		18,370
9	1	17,069	33,458	16,275	0,02	0,238	0,590	0,22	3,581	6,488		
	2	17,072	32,971	15,768	0,00	0,000	0,000	0,295	0,05	0,789	1,430	3,959
	3	17,137	33,008	15,761	0,03	0,473	0,867	0,28	4,413	7,996		
	4	17,065	32,888	15,713	0,04	0,629	1,139	0,996	0,49	7,699	13,951	10,974
	5	17,162	33,284	16,020	5,60	82,916	168,362	0,45	7,209	13,063		
	6	17,095	32,863	15,460	5,88	81,058	164,997	0,46	7,111	12,886		12,974
	7	17,043	32,447	15,297	6,32	96,677	175,178	0,49	7,960	14,426		
	8	17,093	32,769	15,597	6,04	94,205	170,699	0,43	6,707	12,152		13,006
12	1	16,841	33,300	16,245	0,05	0,812	1,472	0,15	2,365	4,266		9,355
	2	16,877	32,756	15,769	0,01	0,158	0,238	0,879	0,15	2,365	4,266	9,355
	3	17,152	32,990	15,718	0,05	0,786	1,424	0,44	6,916	12,532		17,501
	4	17,089	32,895	15,897	0,06	1,256	2,276	0,79	12,401	22,470		17,501
	5	17,066	33,224	16,026	7,03	112,682	204,143	0,60	9,616	17,423		14,028
	6	16,987	32,537	15,442	4,95	76,437	138,505	0,38	5,898	10,633		
	7	16,878	32,247	15,262	8,07	123,168	223,176	0,70	10,864	19,359		
	8	17,293	32,999	15,597	4,81	75,021	135,958	0,39	6,083	11,022		15,190

15	1	17,050	33,415	16,251	0.08	1,300	2,356	1,750	0.58	9,426	17,079
	2	17,230	33,132	15,791	0.04	0.632	1,145	1,750	0.34	5,369	9,729
	3	17,235	33,082	15,737	0.17	2,675	4,948	1,750	1.20	18,884	34,218
	4	17,340	33,167	15,717	0.13	2,043	3,702	4,275	1.04	16,346	29,318
	5	17,254	33,376	16,010	3.65	56,595	103,278	117,274	0.30	4,803	8,703
	6	16,679	32,401	15,414	4.70	72,446	131,273	117,274	0.38	5,649	10,055
	7	17,091	32,476	15,278	5.02	78,896	138,973	117,255	0.43	6,570	11,904
	8	17,301	33,003	15,593	4.15	64,710	117,255	128,114	0.35	5,457	9,889
18	1	16,755	33,089	16,219	0.12	1,946	3,527	2,905	0.89	11,191	20,279
	2	16,708	32,561	15,745	0.08	1,280	2,322	2,905	0.51	6,030	14,550
	3	17,067	32,910	15,733	0.16	2,517	4,361	4,977	0.84	13,216	23,947
	4	16,844	32,417	15,663	0.19	2,976	5,393	4,977	1.20	18,796	34,058
	5	17,139	33,226	15,975	3.94	62,942	114,051	116,885	0.31	4,952	8,974
	6	16,988	32,533	15,437	4.28	69,070	119,719	116,885	0.36	5,557	10,070
	7	16,758	32,122	15,266	4.02	61,370	111,202	107,073	0.35	5,343	9,692
	8	17,198	32,572	15,565	3.65	56,812	102,944	107,073	0.34	5,292	9,589
21	1	17,058	32,627	16,118	0.21	3,385	6,138	4,908	1.08	17,408	31,543
	2	17,175	32,554	15,857	0.13	2,053	3,893	4,908	0.69	10,769	19,550
	3	17,124	32,859	15,626	0.22	3,438	6,229	7,054	1.02	15,938	28,880
	4	16,960	32,564	15,595	0.25	3,899	7,054	6,647	1.31	20,429	37,018
	5	17,110	33,120	15,899	3.60	57,235	103,710	101,384	0.35	5,596	10,083
	6	17,040	32,454	15,307	3.57	54,945	99,018	101,384	0.35	5,507	9,708
	7	17,224	32,489	15,159	3.12	47,298	85,700	89,379	0.33	5,002	9,084
	8	17,244	32,821	15,469	3.32	51,356	93,057	89,379	0.34	5,259	9,530
24	1	17,167	33,397	16,097	0.24	3,863	7,000	6,330	1.03	16,580	30,043
	2	17,224	32,953	15,820	0.20	3,124	5,661	6,330	0.92	14,370	26,039
	3	17,135	32,842	15,598	0.29	4,523	8,196	8,488	1.17	18,249	33,068
	4	17,061	32,729	15,559	0.31	4,823	8,740	8,488	1.40	21,783	39,470
	5	17,168	33,163	15,874	2.90	46,034	83,414	87,333	0.28	4,445	8,054
	6	17,100	32,514	15,307	3.29	50,390	91,251	87,333	0.33	5,051	9,153
	7	16,988	32,244	15,150	2.83	42,874	77,688	78,213	0.29	4,484	8,126
	8	17,085	32,657	15,484	2.81	43,453	78,737	78,213	0.31	4,696	8,510
27	1	16,864	33,063	15,096	0.40	6,485	11,559	9,368	1.54	24,773	44,889
	2	17,205	32,938	15,624	0.25	3,906	7,078	9,368	1.03	16,092	29,159
	3	16,647	32,940	15,584	0.34	5,289	9,584	11,146	1.67	26,025	47,158
	4	17,273	32,939	15,557	0.45	7,013	12,707	11,146	1.40	21,780	39,466
	5	17,220	33,197	16,868	3.23	51,247	92,860	81,673	0.32	5,077	9,200
	6	17,098	32,490	16,315	2.54	38,900	70,488	81,673	0.25	3,829	6,988
	7	17,177	32,434	16,151	3.13	47,422	85,929	75,310	0.35	5,303	9,609
	8	17,055	32,618	15,455	2.31	35,701	64,690	75,310	0.23	3,555	6,441

30	1	17,290	33,506	16,103	0.31	4,992	9,046	1,09	17,553	31,805	33,095	32,450
	2	16,867	32,587	15,611	0.31	4,839	8,769	1,17	18,265	33,095	33,095	32,450
	3	17,150	32,855	15,966	0.42	6,550	11,869	1,34	20,988	37,868	41,482	39,665
	4	16,463	32,138	15,596	0.40	8,226	11,282	1,47	22,852	37,868	41,482	39,665
	5	17,135	33,110	15,854	2.64	41,881	75,888	0.25	3,966	7,186	7,186	7,611
	6	17,142	32,641	15,292	2.94	44,958	81,465	0.29	4,435	8,035	8,035	7,611
	7	16,969	32,204	15,129	2.29	34,546	62,778	0.24	3,681	6,578	6,578	6,511
	8	17,270	32,855	16,457	2.29	35,396	64,138	0.23	3,555	6,442	6,442	6,511
33	1	17,054	33,241	16,074	0.34	5,465	9,903	1,15	16,486	33,496	33,496	33,496
	2	17,030	32,761	15,612	0.35	5,484	9,901	1,21	18,890	34,229	33,863	33,863
	3	17,021	32,720	15,550	0.46	7,015	12,712	1,35	21,046	38,136	40,298	39,217
	4	17,301	32,962	15,562	0.42	6,532	11,636	1,43	22,240	40,298	40,298	39,217
	5	17,077	33,067	15,869	2.48	39,355	71,311	0.22	3,491	6,326	6,326	6,494
	6	17,112	32,537	15,318	2.75	42,124	78,328	0.24	3,678	6,661	6,661	6,494
	7	17,086	32,346	15,154	2.60	42,431	76,885	0.30	4,548	8,238	8,238	8,238
	8	18,981	32,553	15,154	2.10	32,474	55,843	0.20	3,093	5,604	5,604	5,621
36	1	17,130	33,348	16,454	0.28	4,509	8,771	0.84	13,523	24,513	24,513	27,817
	2	17,068	32,821	15,814	0.31	4,840	8,771	1,10	17,175	31,121	31,121	27,817
	3	17,133	32,846	15,604	0.48	7,022	12,723	1,28	19,973	36,191	36,191	37,854
	4	17,012	32,869	15,578	0.40	6,231	11,291	1,40	21,809	39,518	39,518	37,854
	5	17,054	33,035	15,870	2.22	36,231	63,839	0.20	3,174	5,751	5,751	5,082
	6	17,216	32,615	15,293	2.88	40,985	74,295	0.23	3,517	6,373	6,373	5,082
	7	17,092	32,968	15,168	2.31	35,038	63,468	0.24	3,540	6,596	6,596	6,602
	8	17,107	32,694	15,469	2.40	36,816	66,710	0.24	3,712	6,727	6,727	6,602
39	1	17,130	33,310	16,068	0.40	6,427	11,648	1,11	17,855	32,317	32,317	33,972
	2	16,894	32,598	15,605	0.41	6,398	11,593	1,11	17,855	32,317	32,317	33,972
	3	17,161	32,854	15,584	0.51	7,948	14,401	1,26	19,862	35,628	35,628	32,656
	4	17,071	32,719	15,539	0.36	5,439	9,855	1,32	20,571	37,274	37,274	32,656
	5	16,759	32,714	15,844	1.67	26,460	47,845	0.18	2,862	5,168	5,168	4,682
	6	16,921	32,321	15,263	1.61	24,622	44,614	0.15	2,294	4,157	4,157	4,682
	7	16,956	32,213	15,152	1.95	23,546	43,538	0.19	2,878	5,217	5,217	5,411
	8	16,988	32,162	15,465	2.03	24,189	43,538	0.20	3,083	5,605	5,605	5,411
42	1	16,921	33,119	16,095	0.42	6,766	12,242	1,11	17,855	32,353	32,353	33,975
	2	17,110	32,811	15,592	0.43	6,704	12,149	1,26	19,846	35,598	35,598	33,975
	3	17,076	32,772	15,588	0.55	8,573	15,535	1,31	20,420	37,001	37,001	36,085
	4	17,028	32,666	15,561	0.49	7,820	13,808	1,39	21,618	39,168	39,168	36,085
	5	17,059	33,046	15,576	2.05	32,548	58,973	0.18	2,868	5,176	5,176	5,361
	6	17,077	32,484	15,300	2.40	36,720	66,538	0.20	3,080	5,545	5,545	5,361
	7	16,830	32,089	15,153	1.86	29,700	53,816	0.19	2,678	5,217	5,217	5,361
	8	17,124	32,980	15,448	1.93	29,614	54,924	0.19	2,835	5,318	5,318	5,268

45	1	17,231	33,423	18,079	0.58	9,826	18,899	16,079	1.48	23,476	42,538
	2	17,188	32,892	15,566	0.54	8,421	15,259	16,079	1.48	23,080	41,822
	3	17,012	32,710	15,569	0.78	12,169	22,033	17,979	1.79	27,904	50,562
	4	16,846	32,508	15,553	0.90	17,777	24,091	18,082	1.37	21,308	38,610
	5	17,229	33,189	15,563	1.64	25,882	47,088	65,250	0.15	2,377	4,308
	6	17,043	32,472	15,322	2.86	43,820	79,402	65,250	0.25	3,984	7,218
	7	16,917	32,174	15,151	1.86	28,181	51,064	51,862	0.19	2,879	5,763
	8	16,937	32,504	15,459	1.88	29,063	52,681	51,862	0.19	2,937	5,269
48	1	17,121	32,297	16,054	0.56	8,990	16,290	15,868	1.33	21,351	38,689
	2	17,033	32,737	15,585	0.49	7,641	13,846	15,088	1.32	20,571	36,174
	3	16,258	28,951	15,584	0.83	8,818	17,790	16,356	1.35	20,975	38,007
	4	17,334	32,980	15,537	0.53	8,235	14,921	16,356	0.18	2,851	5,166
	5	16,992	32,943	15,840	1.89	29,938	64,247	57,451	0.21	2,904	5,214
	6	17,235	32,627	15,285	2.19	33,474	60,555	57,451	0.21	2,904	5,214
	7	17,271	32,491	15,114	1.89	30,077	54,500	59,587	0.19	2,838	5,319
	8	17,118	32,677	15,451	1.88	29,048	52,034	59,587	0.19	2,838	5,319
51	1	16,971	32,878	15,798	0.62	9,794	17,748	15,258	1.49	23,537	42,648
	2	17,208	32,633	15,320	0.46	7,047	12,769	15,258	1.22	18,690	33,886
	3	16,817	32,348	15,324	0.55	9,960	18,048	15,815	1.36	20,840	37,763
	4	16,955	32,359	15,207	0.48	7,495	13,582	15,815	1.30	19,888	36,033
	5	16,957	32,339	15,573	1.86	25,651	46,842	50,542	0.26	4,049	7,337
	6	16,943	32,091	15,018	1.99	29,938	54,242	50,542	0.30	4,513	8,177
	7	16,703	31,653	14,878	1.74	25,884	46,902	45,468	0.25	3,797	7,817
	8	16,994	32,289	15,189	1.60	24,302	44,035	45,468	0.25	3,797	7,817
54	1	17,029	32,929	15,789	0.65	10,203	18,597	15,823	1.48	23,368	42,344
	2	16,938	32,369	15,327	0.47	7,202	13,049	15,823	1.19	18,234	33,040
	3	17,242	32,678	15,327	0.71	10,882	18,718	20,089	1.87	28,533	51,701
	4	17,006	30,071	15,258	0.74	11,291	20,459	20,089	0.00	0,000	0,000
	5	16,986	32,669	15,558	1.95	30,350	54,994	53,883	0.28	4,217	7,641
	6	17,086	32,231	15,090	1.92	28,914	52,383	53,883	0.26	3,875	7,022
	7	17,187	32,196	14,905	1.55	23,102	41,861	42,708	0.25	3,603	6,967
	8	17,233	32,553	15,214	1.58	24,037	43,556	42,708	1.16	18,331	33,215
57	1	16,902	32,815	15,802	0.51	8,059	14,503	13,419	1.11	17,032	30,863
	2	17,234	32,686	15,345	0.44	6,752	12,234	13,419	1.37	21,017	38,082
	3	17,151	32,599	15,341	0.71	10,892	19,736	16,668	1.24	18,964	34,364
	4	16,850	32,251	15,294	0.49	7,494	13,579	16,668	0.37	5,769	10,453
	5	16,890	32,590	15,581	1.99	31,028	56,219	44,761	0.21	3,164	5,732
	6	16,504	31,874	15,085	1.22	19,379	33,302	44,761	0.19	2,833	5,133
	7	17,102	32,118	14,910	0.99	14,761	26,745	22,994	0.13	1,972	3,574
	8	16,743	32,020	15,171	0.70	10,620	19,243	22,994	0.13	1,972	3,574

60	1	16.961	32.870	15.808	0.27	4.288	7.734	7.082	0.58	9.168	16.614	15.809
	2	17.105	32.547	15.355	0.23	3.927	6.391	7.082	0.54	8.281	15.005	15.809
	3	17.156	32.622	15.359	0.32	4.915	8.906	7.920	0.61	8.755	15.864	16.383
	4	17.103	32.617	15.359	0.25	3.827	6.935	7.920	0.61	9.338	18.921	16.383
	5	17.084	32.787	15.844	0.68	10.804	18.214	20.918	0.12	1.871	3.361	3.603
	6	17.049	32.186	15.042	0.83	12.485	22.622	20.918	0.14	2.108	3.316	3.603
	7	16.993	32.006	14.908	0.68	10.137	18.369	18.411	0.13	1.938	3.512	3.603
	8	17.108	32.415	15.201	0.67	10.184	18.454	18.411	0.13	1.978	3.581	3.603
63	1	17.163	33.081	15.807	0.23	4.426	8.020	7.208	0.57	9.010	16.326	15.807
	2	17.025	32.478	15.346	0.29	3.528	6.395	7.208	0.53	8.133	14.797	15.807
	3	17.136	32.563	15.346	0.32	4.902	8.933	7.907	0.58	8.579	15.545	15.807
	4	17.104	32.511	15.300	0.25	3.825	6.931	7.907	0.53	8.874	16.090	15.812
	5	16.849	32.534	15.578	0.81	12.190	22.089	20.358	0.14	2.107	3.818	3.602
	6	16.974	32.129	15.050	0.81	12.190	22.089	20.358	0.14	2.107	3.818	3.602
	7	16.655	15.72	14.905	0.68	8.638	17.828	17.868	0.12	1.824	3.306	3.408
	8	17.165	32.473	15.204	0.65	8.682	17.907	17.868	0.12	1.824	3.306	3.408
66	1	17.088	32.984	15.705	0.29	4.160	9.172	9.057	0.55	8.440	15.293	18.666
	2	16.922	32.375	15.348	0.23	3.836	6.952	9.057	0.53	8.130	14.732	18.666
	3	17.115	32.562	15.348	0.25	4.160	9.172	9.057	0.53	8.130	14.732	18.666
	4	16.903	32.021	15.311	0.33	3.602	6.952	9.057	0.53	8.130	14.732	18.666
	5	17.196	32.888	15.583	1.29	20.102	36.425	39.583	0.23	3.584	6.494	21.854
	6	17.124	32.253	15.824	1.57	23.597	42.740	39.583	0.27	4.056	7.950	8.922
	7	16.883	31.869	15.209	1.23	18.305	33.168	33.532	0.24	3.572	6.494	8.922
	8	17.174	32.459	15.209	1.23	18.305	33.168	33.532	0.24	3.572	6.494	8.922
69	1	16.219	32.895	15.805	0.55	8.693	15.752	14.404	1.04	18.438	29.785	8.405
	2	17.089	32.927	15.331	0.47	7.205	13.056	14.404	1.01	15.484	28.057	28.821
	3	17.259	32.729	15.362	0.73	11.216	20.331	17.091	1.16	17.820	32.291	28.821
	4	16.957	32.264	15.300	0.50	7.650	13.952	17.091	1.11	16.983	30.773	31.532
	5	17.196	32.863	15.558	1.28	18.603	35.521	37.968	0.23	3.578	6.484	8.405
	6	17.075	32.272	15.091	1.48	22.385	40.471	37.968	0.26	3.880	7.030	8.405
	7	17.161	32.188	14.923	1.36	20.145	36.504	34.508	0.22	3.945	6.081	6.548
	8	17.244	32.554	15.204	1.18	17.940	32.509	34.508	0.22	3.945	6.081	6.548
72	1	16.950	32.857	15.798	0.60	9.478	17.174	15.117	1.10	17.376	31.465	29.498
	2	16.787	32.230	15.336	0.47	7.208	13.060	15.117	0.99	15.182	27.510	29.498
	3	17.186	32.829	15.336	0.86	10.122	18.340	15.831	1.06	18.237	29.421	29.498
	4	17.196	32.821	15.318	0.48	7.353	13.323	15.831	1.06	18.237	29.421	29.498
	5	17.192	32.892	15.591	1.18	18.387	33.338	35.653	0.21	3.274	5.963	6.244
	6	17.092	32.273	15.075	1.39	20.955	37.970	35.653	0.24	3.618	6.568	6.244
	7	17.168	32.166	14.914	1.18	17.598	31.988	31.634	0.22	3.261	5.946	6.244
	8	16.956	32.256	15.192	1.14	17.318	31.381	31.634	0.21	3.190	5.761	5.863